

DICOM Conformance Statement
for
NW500

Ver.2
Oct. 26th, 2022
TOPCON Corporation
ALL RIGHTS RESERVED

1. CONFORMANCE STATEMENT OVERVIEW

This document declares conformance to DICOM V3.0 for NW500. NW500 works as an Acquisition Modality and allows acquisition and storage of acquisition data. The following table indicates DICOM SOP Classes which NW500 supports.

SOP Classes	User of Service (SCU)	Provider of Service (SCP)
VERIFICATION		
Verification	Yes	No
TRANSFER		
Ophthalmic Photography 8bit Image Storage	Yes	No
VL Photographic Image Storage	Yes	No
Secondary Capture Image Storage	Yes	No
WORKFLOW MANAGEMENT		
Modality Worklist Information Model - FIND	Yes	No
Modality Performed Procedure Step SOP Class	Yes	No
Storage Commitment Push Model SOP Class	Yes	No
QUERY / RETRIEVE		
Patient Root Query / Retrieve Information Model – FIND	Yes	No

2. TABLE OF CONTENTS

1. CONFORMANCE STATEMENT OVERVIEW	2
2. TABLE OF CONTENTS	3
3. INTRODUCTION	4
3.1. REVISION HISTORY	4
3.2. AUDIENCE	4
3.3. ABBREVIATIONS	4
3.4. REFERENCES	4
4. NETWORKING	5
4.1. IMPLEMENTATION MODEL	5
4.1.1. Application Data Flow	5
4.1.2. Function Definitions of NW500 AE	6
4.2. NW500 AE SPECIFICATIONS	7
4.2.1. Verification	7
4.2.2. Worklist	9
4.2.3. Storage	12
4.2.4. Storage Commitment	14
4.2.5. Patient Root Query	16
4.2.1. MPPS	18
4.3. NETWORK INTERFACES	23
4.3.1. Physical Network Interface	23
4.3.2. IPv4 and IPv6 Support	23
4.4. DATA DICTIONARY OF PRIVATE ATTRIBUTES	23
4.5. STANDARD EXTENDED / SPECIALIZED / PRIVATE SOP CLASSES	23
4.6. PRIVATE TRANSFER SYNTAXES	23
4.7. CONFIGURATION	24
4.7.1. AE Title / Presentation Address Mapping	24
4.7.2. Parameters	24
5. MEDIA INTERCHANGE	25
6. SUPPORT OF CHARACTER SETS	25
6.1. Without Code Extensions	25
6.2. With Code Extensions	25
ANNEX A	26

3. INTRODUCTION

3.1. REVISION HISTORY

Document Version	Date of Issue	Author	Description
1.00	Jun. 7th, 2022	Eye Care Product System Dept.	The Publication of the first edition.
2.00	Oct. 26th, 2022	Eye Care Product System Dept.	<ul style="list-style-type: none">• Added MPPS function.• Deleted unnecessary definitions.• Added the definitions to use to store images.

3.2. AUDIENCE

This document is intended for hospital staffs, health system integrators, software engineers, service staffs who have a basic knowledge of DICOM.

3.3. ABBREVIATIONS

- ACR American College of Radiology
- AE Application Entity
- ANSI American National Standards Institute
- ASCII American Standard Code for Information Interchange
- DICOM Digital Imaging and Communication in Medicine
- DIMSE DICOM Message Service Element
- DIMSE-C DICOM Message Service Element-Composite
- DIMSE-N DICOM Message Service Element-Normalized
- IE Information Entity
- IOD Information Object Definition
- ISO International Standards Organization
- NEMA National Electrical Manufacture Association
- OSI Open Systems Interconnection
- PDU Protocol Data Unit
- SCP Service Class Provider
- SCU Service Class User
- SOP Service Object-Pair
- TCP/IP Transmission Control Protocol/Internet Protocol
- UID Unique Identifier


3.4. REFERENCES

- ACR-NEMA DICOM (Digital Imaging and Communication in Medicine)

4. NETWORKING

4.1. IMPLEMENTATION MODEL

4.1.1. Application Data Flow



4.1.2. Function Definitions of NW500 AE

4.1.2.1. Verification Function Definitions

NW500 AE checks the connection settings of DICOM by sending a C-ECHO request message to a Verification AE.

4.1.2.2. Worklist Function Definitions

NW500 AE receives a study order list by sending a search condition with a C-FIND request message to a remote Worklist AE. And the study order list is displayed on NW500.

4.1.2.3. Storage Function Definitions

NW500 AE sends a captured image with a C-STORE request message to a remote Storage AE. NW500 AE corresponds to OP, VL and SC modality.

4.1.2.4. Storage Commitment Function Definitions

NW500 AE sends a storage commitment message with a N-ACTION request message to a remote Storage Commitment AE. After a remote Storage Commitment AE completes the requested process, a NW500 AE receives the result of the storage commitment with a N-EVENT-REPORT message from a remote Storage Commitment AE.

4.1.2.5. Patient Root Query Function Definitions

NW500 AE receives a patient list by sending a condition of query with a C-FIND request message to a remote Patient Root Query AE. And the patient list is displayed on NW500.

4.1.2.6. MPPS Function Definitions

NW500 AE sends a status of study implementation with a N-CREATE request message or a N-SET request message to a remote MPPS AE.

4.2. NW500 AE SPECIFICATIONS

4.2.1. Verification

NW500 AE provides Standard Conformance to the following DICOM V3.0 SOP class as a SCU.

4.2.1.1. SOP Class

SOP Class Name	SOP Class UID	Role
Verification	1.2.840.10008.1.1	SCU

4.2.1.2. Association Policies

4.2.1.2.1. General

The Application Context Name for DICOM 3.0 is the only Application Context proposed.

Application Context Name	1.2.840.10008.3.1.1.1
--------------------------	-----------------------

4.2.1.2.2. Number of Associations

NW500 AE can establish only one association simultaneously.

4.2.1.2.3. Asynchronous Nature

NW500 AE allows only a single operation for an association. Therefore, an asynchronous operation is not supported.

4.2.1.2.4. Implementation Identifying Information

NW500 AE specifies the following implementation identifying information.

Implementation Class UID	1.2.392.200106.1650.1.5
Implementation Version Name	TP_VER_NW500_100

4.2.1.2.5. Association Initiation Policy

4.2.1.2.5.1. Description and Sequencing Activities

After NW500 AE establishes a new association to check the connection settings, it sends a verification request message to a remote Verification AE by using a C-ECHO service of DIMSE-C.

4.2.1.2.5.2. Proposed Presentation Contexts

NW500 AE can propose the Presentation Contexts shown in the following table.

Presentation Context					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID List		
Verification	1.2.840.10008.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None

4.2.1.2.5.3. SOP Specific Conformance for Verification SOP Class

When the result of data communications is successful, NW500 AE judges that verification is successful even if the release of an association fails.

4.2.1.2.6. A receiving policy of association

NW500 AE does not accept any association requests.

4.2.2. Worklist

NW500 AE provides Standard Conformance to the following DICOM V3.0 SOP class as a SCU.

4.2.2.1. SOP Class

SOP Class Name	SOP Class UID	Role
Modality Worklist Information Model-FIND	1.2.840.10008.5.1.4.31	SCU

4.2.2.2. Association Policies

4.2.2.2.1. General

The Application Context Name for DICOM 3.0 is the only Application Context proposed.

Application Context Name	1.2.840.10008.3.1.1.1

4.2.2.2.2. Number of Associations

NW500 AE can establish only one association simultaneously.

4.2.2.2.3. Asynchronous Nature

NW500 AE allows only a single operation for an association. Therefore, an asynchronous operation is not supported.

4.2.2.2.4. Implementation Identifying Information

NW500 AE specifies the following implementation identifying information.

Implementation Class UID	1.2.392.200106.1650.1.1
Implementation Version Name	TP_MWL_NW500_100

4.2.2.2.5. Activity

4.2.2.2.5.1. Description and Sequencing Activities

After NW500 AE establishes a new association to receive a study order list, it sends a search condition to a remote Worklist AE by using a C-FIND service of DIMSE-C.

4.2.2.2.5.2. Proposed Presentation Contexts

NW500 AE can propose the Presentation Contexts shown in the following table.

Presentation Context					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID List		
Modality Worklist Information Model-FIND	1.2.840.10008.5.1.4.31	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None

4.2.2.2.5.3. SOP Specific Conformance for Worklist SOP Class

When the result of data communications is successful, NW500 AE judges that the acquisition of a study order list is successful even if the release of an association fails.

Tag Description	Tag	Matching	Display
Scheduled Procedure Step			
Scheduled Procedure Step Sequence	(0040,0100)		
>Scheduled Station AE Title	(0040,0001)	S+	
>Scheduled Procedure Step Start Date	(0040,0002)	S	
>Scheduled Procedure Step Start Time	(0040,0003)		
>Modality	(0008,0060)	S+	
>Scheduled Performing Physicians Name	(0040,0006)		
>Scheduled Procedure Step Description	(0040,0007)		
>Scheduled Station Name	(0040,0010)		
>Scheduled Procedure Step Location	(0040,0011)		
>Scheduled Action Item Code Sequence	(0040,0008)		
>>Code Value	(0008,0100)		
>>Coding Scheme Designator	(0008,0102)		
>>Coding Scheme Version	(0008,0103)		
>>Coding Meaning	(0008,0104)		
>Pre-Medication	(0040,0012)		
>Scheduled Procedure Step ID	(0040,0009)		
>Requested Contrast Agent	(0032,1070)		
>Scheduled Procedure Step Status	(0040,0020)		
Requested Procedure			
Requested Procedure ID	(0040,1001)		
Requested Procedure Description	(0032,1060)		
Requested Procedure Code Sequence	(0032,1064)		
Study Instance UID	(0020,000D)		
Study Date	(0008,0020)		
Study Time	(0008,0030)		
Requested Procedure Priority	(0040,1003)		
Patient Transport Arrangements	(0040,1004)		
Imaging Service Request			
Accession Number	(0008,0050)		○
Requesting Physician	(0032,1032)		
Referring Physician's Name	(0008,0090)		
Requesting Service	(0032,1033)		
Visit Identification			
Admission ID	(0038,0010)		
Visit Status			
Current Patient Location	(0038,0300)		
Patient Identification			
Patient's Name	(0010,0010)	*	○
Patient ID	(0010,0020)	*	○
Patient Demographic			
Patients Birth Date	(0010,0030)		
Patient's Sex	(0010,0040)		○
Patient's Weight	(0010,1030)		
Ethnic Group	(0010,2160)		
Confidentiality Constraint On Patient Data	(0040,3001)		
Patient State	(0038,0500)		
Pregnancy Status	(0010,21C0)		
Medical Alerts	(0010,2000)		
Contrast Allergies	(0010,2110)		
Special Needs	(0038,0050)		
Referenced Study Sequence	(0008,1110)		

Tag : Tag Number

Matching : Search key for updating the worklist. 'S' provides an attribute for a single inspection inspection. '+' indicates a configurable item in the setting page.
*' indicates a wildcard search is available.

Display : '×' indicates the items that NW500 can display on the monitor screen.

4.2.2.2.6. A receiving policy of association

NW500 AE does not accept any association requests.

4.2.3. Storage

NW500 AE provides Standard Conformance to the following DICOM V3.0 SOP class as a SCU.

4.2.3.1. SOP Classes

SOP Class Name	SOP Class UID	Role
Ophthalmic Photography 8 Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.1	SCU
VL Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4	SCU
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	SCU

4.2.3.2. Association Policies

4.2.3.2.1. General

The Application Context Name for DICOM 3.0 is the only Application Context proposed.

Application Context Name	1.2.840.10008.3.1.1.1

4.2.3.2.2. Number of Associations

NW500 AE can establish only one association simultaneously.

4.2.3.2.3. Asynchronous Nature

NW500 AE allows only a single operation for an association. Therefore, an asynchronous operation is not supported.

4.2.3.2.4. Implementation Identifying Information

NW500 AE specifies the following implementation identifying information.

Implementation Class UID	1.2.392.200106.1650.1.2
Implementation Version Name	TP_STO_NW500_100

4.2.3.2.5. Activity

4.2.3.2.5.1. Description and Sequencing Activities

After NW500 AE establishes a new association to store an image, it sends an image to a remote Storage AE by using a C-STORE service of DIMSE-C.

4.2.3.2.5.2. Proposed Presentation Contexts

The store AE can propose the Presentation Contexts shown in the following table.

Presentation Context					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
VL Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCU	None

		JPEG Baseline Lossy Compression (*)	1.2.840.10008.1.2.4.50	SCU	None
SC Image Storage	1.2.840.10008.5.1.4.1.1.7	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCU	None
		JPEG Baseline Lossy Compression (*)	1.2.840.10008.1.2.4.50	SCU	None
Ophthalmic Photography 8 Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCU	None
		JPEG Baseline Lossy Compression (*)	1.2.840.10008.1.2.4.50	SCU	None

(*) JPEG Baseline (Process 1)

4.2.3.2.5.3. SOP Specific Conformance for Storage SOP Class

When the result of data communications is successful, NW500 AE judges that the storing of an image is successful even if the release of an association fails.

The details of IODs which is sent by NW500 AE are described in ANNEX A.

4.2.3.2.6. A receiving policy of association

The store AE does not accept any association requests.

4.2.4. Storage Commitment

NW500 AE provides Standard Conformance to the following DICOM V3.0 SOP class as a SCU.

4.2.4.1. SOP Class

SOP Class Name	SOP Class UID	Role
Storage Commitment Push Model SOP Class	1.2.840.10008.1.20.1	SCU

4.2.4.2. Association Policies

4.2.4.2.1. General

The Application Context Name for DICOM 3.0 is the only Application Context proposed.

Application Context Name	UID
	1.2.840.10008.3.1.1.1

4.2.4.2.2. Number of Associations

NW500 AE can establish only one association simultaneously.

4.2.4.2.3. Asynchronous Nature

NW500 AE allows only a single operation for an association. Therefore, an asynchronous operation is not supported.

4.2.4.2.4. Implementation Identifying Information

NW500 AE specifies the following implementation identifying information.

Implementation Class UID	1.2.392.200106.1650.1.4
Implementation Version Name	TP_COM_NW500_100

4.2.4.2.5. Activity

4.2.4.2.5.1. Description and Sequencing Activities

After NW500 AE establishes a new association to ask to save an image enduringly, it sends a storage commitment request message to a remote Storage Commitment AE by using a N-ACTION service of DIMSE-N. Then, when a remote Storage Commitment AE completes a process of a storage commitment request, it sends a result of storage commitment to NW500 AE by using a N-EVENT-REPORT service of DICOM-N. Therefore, NW500 AE accepts a transmission request of a result in the association which was established for N-ACTION. However, if the association already is released at that time, a remote Storage Commitment AE establishes a new association to send a transmission request of a result.

4.2.4.2.5.2. Proposed Presentation Contexts

NW500 AE can propose the Presentation Contexts shown in the following table.

Presentation Context					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID List		
Storage Commitment Push Model SOP Class	1.2.840.10008.1.20.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None

4.2.4.2.5.3. SOP Specific Conformance for Storage Commitment SOP Class

When the result of data communications is successful, NW500 AE judges that a request of saving an image enduringly is successful even if the release of an association fails.

Attribute Name	Tag	Type	Attribute Description
Storage Commitment Request			
Transaction UID	(0008,1195)	1	A value created by NW500.
Referenced SOP Sequence	(0008,1199)	1	A value acquired from Storage data.
>Referenced SOP Class UID	(0008,1150)	1	A value acquired from Storage data.
>Referenced SOP Instance UID	(0008,1155)	1	A value acquired from Storage data.

4.2.4.2.6. A receiving policy of association

NW500 AE accepts only the following association request which is requested from a remote Storage Commitment AE.

- A notification request of a result of a storage commitment

4.2.5. Patient Root Query

NW500 AE provides Conformance to the following DICOM V3.0 SOP class as a SCU.

4.2.5.1. SOP Class

SOP Class Name	SOP Class UID	Role
Patient Root Query / Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.2.1	SCU

4.2.5.2. Association Policies

4.2.5.2.1. General

The Application Context Name for DICOM 3.0 is the only Application Context proposed.

Application Context Name	1.2.840.10008.3.1.1.1

4.2.5.2.2. Number of Associations

NW500 AE can establish only one association simultaneously.

4.2.5.2.3. Asynchronous Nature

NW500 AE allows only a single operation for an association. Therefore, an asynchronous operation is not supported.

4.2.5.2.4. Implementation Identifying Information

Implementation Class UID	1.2.392.200106.1650.1.6
Implementation Version Name	TP_QUE_NW500_100

4.2.5.2.5. Activity

4.2.5.2.5.1. Description and Sequencing Activities

After NW500 AE establishes a new association to receive a patient information list from PACS, it sends a search condition to a remote Patient Root Query AE by using a C-FIND service of DIMSE-C.

4.2.5.2.5.2. Proposed Presentation Contexts

NW500 AE can propose the Presentation Contexts shown in the following table.

Presentation Context					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID List		
Patient Root Query / Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.2.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None

4.2.5.2.5.3. SOP Specific Conformance for Patient Root Query SOP Class

When the result of data communications is successful, NW500 AE judges that the acquisition of a patient information list is successful even if the release of an association fails.

And the value of Query / Retrieve level is fixed to 'PATIENT'.

Tag Description	Tag	Matching	Display
Patient Level			
Patient's Name	(0010,0010)		×
Patient ID	(0010,0020)	*	×
Patients Birth Date	(0010,0030)		
Patient's Sex	(0010,0040)		×
Ethnic Group	(0010,2160)		

Tag : Tag Number

Matching : Search key for updating the patient list. '*' indicates a wildcard search is available.

Display : '×' indicates the items that NW500 can display on the monitor screen.

4.2.5.2.6. A receiving policy of association

NW500 AE does not accept any association requests.

4.2.1. MPPS

NW500 AE provides Standard Conformance to the following DICOM V3.0 SOP class as a SCU.

4.2.1.1. SOP Classes

SOP Class Name	SOP Class UID	Role
Modality Performed Procedure Step	1.2.840.10008.3.1.2.3.3	SCU

4.2.1.2. Association Policies

4.2.1.2.1. General

The Application Context Name for DICOM 3.0 is the only Application Context proposed.

DICOM V3.0 Application Context Name	UID
	1.2.840.10008.3.1.1.1

4.2.1.2.2. Number of Associations

NW500 AE can establish only one association simultaneously.

4.2.1.2.3. Asynchronous

NW500 AE allows only a single operation for an association. Therefore, an asynchronous operation is not supported.

4.2.1.2.4. Implementation Identifying Information

NW500 AE specifies the following implementation identifying information.

Implementation Class UID	1.2.392.200106.1650.1.3
Implementation Version Name	TP_MPP_NW500_100

4.2.1.2.5. Activity

4.2.1.2.5.1. Description and Sequencing Activities

After NW500 AE establishes a new association, it sends MPPS information to a remote MPPS AE by using a N-CREATE service or a N-SET service of DIMSE-N to notify status information of “Start study”, “Discontinue study”, or “Complete study”.

4.2.1.2.5.2. Proposed Presentation Contexts

NW500 AE can propose the Presentation Contexts shown in the following table.

Presentation Context					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID		
Modality Performed Procedure Step	1.2.840.10008.3.1.2.3.3	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None

4.2.1.2.5.3. SOP Specific Conformance for MPPS

When the result of data communications is successful, NW500 AE judges that the sending of a status notification is successful even if the release of an association fails.

4.2.1.2.5.3.1. N-CREATE

Attribute Name	Tag	Req. Type N-CREATE	Value	
			With MWL	Without MWL
Specific Character Set	(0008,0005)	1C (Required if an extended or replacement character set is used)	An acquired value from MWL.	No tag / An acquired value by Patient Root Query
Performed Procedure Step Relationship				
Scheduled Step Attribute Sequence	(0040,0270)	1	---	---
>Study Instance UID	(0020,000D)	1	An acquired value from MWL.	A generated value by NW500.
>Referenced Study Sequence	(0008,1110)	2	Empty	Empty
>Accession Number	(0008,0050)	2	An acquired value from MWL.	Empty
>Placer Order Number/Imaging Service Request	(0040,2016)	3	An acquired value from MWL. / No tag	No tag
>Filler Order Number/Imaging Service Request	(0040,2017)	3	An acquired value from MWL. / No tag	No tag
>Requested Procedure ID	(0040,1001)	2	An acquired value from MWL. / Empty	Empty
>Requested Procedure Description	(0032,1060)	2	An acquired value from MWL. / Empty	Empty
>Scheduled Procedure Step ID	(0040,0009)	2	An acquired value from MWL. / Empty	Empty
>Scheduled Procedure Step Description	(0040,0007)	2	An acquired value from MWL. / Empty	Empty
>Scheduled Protocol Code Sequence	(0040,0008)	2	Empty	Empty
Patient's Name	(0010,0010)	2	An acquired value from MWL.	An input value / Empty / An acquired value by Patient Root Query
Patient ID	(0010,0020)	2	An acquired value from MWL.	An input value / Empty
Patient's Birth Date	(0010,0030)	2	An acquired value from MWL.	Empty
Patient's Sex	(0010,0040)	2	An acquired value from MWL.	Empty
Referenced Patient Sequence	(0008,1120)	2	An acquired value from MWL. / Empty	Empty
>Referenced SOP Class UID	(0008,1150)	1	An acquired value from MWL. / No tag	No tag
>Referenced Instance UID	(0008,1155)	1	An acquired value from MWL. / No tag	No tag
Performed Procedure Step Information				

Performed Procedure Step ID	(0040,0253)	1	A generated value by NW500.	A generated value by NW500.
Performed Station AE Title	(0040,0241)	1	A preset value	A preset value
Performed Station Name	(0040,0242)	2	A preset value / Empty	A preset value / Empty
Performed Location	(0040,0243)	2	Empty	Empty
Performed Procedure Step Start Date	(0040,0244)	1	Study start date	Study start date
Performed Procedure Step Start Time	(0040,0245)	1	Study start time	Study start time
Performed Procedure Step Status	(0040,0252)	1	IN PROGRESS	IN PROGRESS
Performed Procedure Step Description	(0040,0254)	2	Empty	Empty
Performed Procedure Type Description	(0040,0255)	2	Empty	Empty
Procedure Code Sequence	(0008,1032)	2	Empty	Empty
Performed Procedure Step End Date	(0040,0250)	2	Empty	Empty
Performed Procedure Step End Time	(0040,0251)	2	Empty	Empty
Image Acquisition Results				
Modality	(0008,0060)	1	A preset value	A preset value
Study ID	(0020,0010)	2	An acquired value from MWL. / A generated value by NW500.	A generated value by NW500.
Performed Protocol Code Sequence	(0040,0260)	2	Empty	Empty
>Code Value	(0008,0100)	1C Shall be present if the code value length is 16 characters or less, and the code value is not a URN or URL.	An acquired value from MWL. / No tag	No tag
>Coding Scheme Designator	(0008,0102)	1C Shall be present if Code Value (0008,0100) or Long Code Value (0008,0119) is present. May be present otherwise.	An acquired value from MWL. / No tag	No tag
>Code Meaning	(0008,0104)	3	An acquired value from MWL. / No tag	No tag
Performed Series Sequence	(0040,0340)	2	---	---
>Performing Physician's Name	(0008,1050)	2	Empty	Empty

>Protocol Name	(0018,1030)	1	NW500 Protocol	NW500 Protocol
>Operators' Name	(0008,1070)	2	Empty	Empty
>Series Instance UID	(0020,000E)	1	A generated value by NW500.	A generated value by NW500.
>Series Description	(0008,103E)	2	Color	Color
>Retrieve AE Title	(0008,0054)	2	Empty	Empty
>Referenced Image Sequence	(0008,1140)	2	Empty	Empty
>Referenced Non-Image Composite SOP Instance Sequence	(0040,0220)	2	Empty	Empty

4.2.1.2.5.3.2. N-SET

Attribute Name	Tag	Req. Type N-SET	Value			
			With MWL		Without MWL	
			Complete	Discontinued	Complete	Discontinued
Specific Character Set	(0008,0005)	1C (Required if an extended or replacement character set is used in an Attribute that is set)	An acquired value from MWL.	An acquired value from MWL.	No tag	No tag
Performed Procedure Step Information						
Performed Procedure Step Status	(0040,0252)	3	COMPLETED	DISCONTINUED	COMPLETED	DISCONTINUED
Performed Procedure Step End Date	(0040,0250)	3	Study end date	No tag	Study end date	No tag
Performed Procedure Step End Time	(0040,0251)	3	Study end time	No tag	Study end time	No tag
Image Acquisition Results						
Performed Protocol Code Sequence	(0040,0260)	3	--- / --- / No tag	No tag	--- / No tag	No tag
>Code Value	(0008,0100)	3 Shall be present if the code value length is 16 characters or less, and the code value is not a URN or URL.	An acquired value from MWL. / An input value / No tag	No tag	An input value / No tag	No tag
>Coding Scheme Designator	(0008,0102)	1C Shall be present if Code Value (0008,0100) or Long Code Value (0008,0119) is present. May be present otherwise.	An acquired value from MWL. / An input value / No tag	No tag	An input value / No tag	No tag
>Code Meaning	(0008,0104)	3	An acquired value from MWL. / An input value	No tag	An input value / No tag	No tag

			/ No tag			
Performed Series Sequence	(0040,0340)	3	---	---	---	---
>Performing Physician's Name	(0008,1050)	2	An acquired value from MWL. / Empty	An acquired value from MWL. / Empty	Empty	Empty
>Protocol Name	(0018,1030)	1	NW500 Protocol	NW500 Protocol	NW500 Protocol	NW500 Protocol
>Operators' Name	(0008,1070)	2	Empty	Empty	Empty	Empty
>Series Instance UID	(0020,000E)	1	A generated value	A generated value	A generated value	A generated value
>Series Description	(0008,103E)	2	color	Color	color	color
>Retrieve AE Title	(0008,0054)	2	Empty	Empty	Empty	Empty
>Referenced Image Sequence	(0008,1140)	2	Find (*1).	Empty	Find (*1).	Empty
>>Referenced SOP Class UID	(0008,1150)	1	SOP Class UID used in the Storage data.	---	SOP Class UID used in the Storage data.	---
>>Referenced SOP Instance UID	(0008,1155)	1	SOP Instance UID used in the Storage data.	---	SOP Instance UID used in the Storage data.	---
>Referenced Non-Image Composite SOP Instance Sequence	(0040,0220)	2	Empty	Empty	Empty	Empty

(*1). Add (0008,1150) and (0008,1155) for the number of DICOM objects stored in one session (from MPSS INPROGRESS to MPSS COMPLETED) using the ITEM tag (FFFE,E000).

4.2.1.2.5.3.3. A receiving policy of association

NW500 AE does not accept any association requests.

4.3. NETWORK INTERFACES

NW500 supports DICOM TCP/IP network communication which is defined in PS 3.8 of DICOM Standard. In addition, NW500 supports the TCP/IP protocol stack of Linux system.

4.3.1. Physical Network Interface

NW500 supports a single network interface. One of the following physical network interfaces will be available.

Ethernet 1000base T
Ethernet 100base T

4.3.2. IPv4 and IPv6 Support

NW500 only supports IPv4 connections.

4.4. DATA DICTIONARY OF PRIVATE ATTRIBUTES

NW500 does not support any Private Attributes.

4.5. STANDARD EXTENDED / SPECIALIZED / PRIVATE SOP CLASSES

NW500 does not support any Extended, Specialized or Private SOP Classes.

4.6. PRIVATE TRANSFER SYNTAXES

NW500 does not support any Private Transfer Syntaxes.

4.7. CONFIGURATION

4.7.1. AE Title / Presentation Address Mapping

NW500 uses the AE titles and TCP/IP ports which are specified in the setting page.

4.7.2. Parameters

Many parameters for the general operations can be specified by using a configuration user interface. The following table shows the configurable parameters for DICOM communication.

SCU		
Items	Parameters	Descriptions
NW500	IP Address	Default: None
	Port Number	Default: 64001
	AE Title	Default: TP_AM_NW500_001
SCP		
Items	Parameters	Descriptions
Modality Worklist	IP Address	Default: None
	Port Number	Default: 0
	AE Title	Default: None
Storage	IP Address	Default: None
	Port Number	Default: 0
	AE Title	Default: None
	Modality	Default: SC The following values are available. - OP - VL - SC
	SC Modality	Default: SC The following values are available. - SC - OP - XC - OT
Transfer Syntax	Default: Implicit VR Little Endian The following values are available. - Implicit VR Little Endian, - Explicit VR Little Endian, - JPEG Baseline Lossy Compression	
Storage Commitment	IP Address	Default: None
	Port Number	Default: 0
	AE Title	Default: None
Patient Root Query	IP Address	Default: None
	Port Number	Default: 0
	AE Title	Default: None
MPPS	IP Address	Default: None
	Port Number	Default: 0
	AE Title	Default: None

5. MEDIA INTERCHANGE

NW500 does not support Media Interchange.

6. SUPPORT OF CHARACTER SETS

NW500 supports the following character sets.

6.1. Without Code Extensions

Defined Term	Character Set
ISO_IR 6	ISO 646
ISO_IR 192	UTF-8

6.2. With Code Extensions

Defined Term	Character Set
ISO 2022 IR 6	ISO 646
ISO 2022 IR 87	JIS X 0208
ISO 2022 IR 13	JIS X 0201

ANNEX A

A.1 IODs of NW500 AE

When NW500 AE sends an image, after generating an IOD which conforms to one of the following SOP classes, NW500 AE sends it.

A.1.1 Ophthalmic Photography 8 Bit Image IOD

Information Entity	Module	Reference	Usage* ¹
Patient	Patient	A.2.1	M
Study	General Study	A.2.2	M
Series	General Series	A.2.3	M
	Ophthalmic Photography Series	A.2.4	M
Frame Of Reference	Synchronization	A.2.5	M
Equipment	General Equipment	A.2.6	M
Image	General Image	A.2.7	M
	Image Pixel	A.2.8	M
	Cine	A.2.14	C
	Multi Frame	A.2.15	M
	Acquisition Context	A.2.9	U
	Ophthalmic Photography Image	A.2.16	M
	Ocular Region Imaged	A.2.17	M
	Ophthalmic Photography Acquisition Parameters	A.2.18	M
	Ophthalmic Photographic Parameters	A.2.19	M
SOP Common	A.2.10	M	

*1: M=Mandatory, C=Conditional, U=User option

A.1.2 VL Photographic Image IOD

Information Entity	Module	Reference	Usage* ¹
Patient	Patient	A.2.1	M
Study	General Study	A.2.2	M
Series	General Series	A.2.3	M
Equipment	General Equipment	A.2.6	M
Image	General Image	A.2.7	M
	Image Pixel	A.2.8	M
	Acquisition Context	A.2.9	M
	VL Image	A.2.12	M
	SOP Common	A.2.10	M

*1: M=Mandatory, C=Conditional, U=User option

A.1.3 Secondary Capture Image IOD

Information Entity	Module	Reference	Usage* ¹
Patient	Patient	A.2.1	M
Study	General Study	A.2.2	M
Series	General Series	A.2.3	M
Equipment	General Equipment	A.2.6	U
	SC Equipment	A.2.11	M
Image	General Image	A.2.7	M
	Image Pixel	A.2.8	M
	SC Image	A.2.13	M
	SOP Common	A.2.10	M

*1: M=Mandatory, C=Conditional, U=User option

A.2 Module lists of IODs

A.2.1 Patient Module

Attribute Name	Tag	Type	Attribute Description
Patient's Name	(0010,0010)	2	A value acquired from MWL, a value acquired from PACS, a value generated by NW500 or empty.
Patient ID	(0010,0020)	2	A value acquired from MWL, a value acquired from PACS or a value generated by NW500.
Patient's Birth Date	(0010,0030)	2	A value acquired from MWL or a value acquired from PACS.
Patient's Sex	(0010,0040)	2	A value acquired from MWL or a value acquired from PACS.

A.2.2 General Study Module

Attribute Name	Tag	Type	Attribute Description
Study Instance UID	(0020,000D)	1	A value acquired from MWL or a value generated by NW500.
Study Date	(0008,0020)	2	A value generated by NW500.
Study Time	(0008,0030)	2	A value generated by NW500 or empty.
Referring Physician's Name	(0008,0090)	2	A value acquired from MWL or empty.
Study ID	(0020,0010)	2	A value acquired from MWL or empty.
Accession Number	(0008,0050)	2	A value acquired from MWL or empty.
Study Description	(0008,1030)	3	fundus photo

A.2.3 General Series Module

Attribute Name	Tag	Type	Attribute Description
Modality	(0008,0060)	1	VL or SC
Series Instance UID	(0020,000E)	1	A value generated by NW500.
Series Number	(0020,0011)	2	A value generated by NW500.
Laterality	(0020,0060)	2C	R or L. In case of OP, this tag does not exist in an IOD.
Series Date	(0008,0021)	3	A value generated by NW500.
Series Time	(0008,0031)	3	A value generated by NW500.
Series Description	(0008,103E)	3	color
Patient Position	(0018,5100)	2C	Empty. In case of VL, this tag exists in an IOD.
Request Attributes Sequence	(0040,0275)	3	-
>Requested Procedure ID	(0040,1001)	1C	A value acquired from MWL.
>Scheduled Procedure Step ID	(0040,0009)	1C	A value acquired from MWL.
> Scheduled Procedure Step Description	(0040,0007)	3	A value acquired from MWL.
> Scheduled Protocol Code Sequence	(0040,0008)	3	-
>>Code Value	(0008,0100)	2	A value acquired from MWL.
>>Coding Scheme Designator	(0008,0102)	2	A value acquired from MWL.
>>Code Meaning	(0008,0104)	2	A value acquired from MWL.

A.2.4 Ophthalmic Photography Series Module

In case of 'Ophthalmic Photography 8 Bit Image', this module is used.

Attribute Name	Tag	Type	Attribute Description
Modality	(0008,0060)	1	OP

A.2.5 Synchronization Module

In case of 'Ophthalmic Photography 8 Bit Image', this module is used.

Attribute Name	Tag	Type	Attribute Description
Synchronization Frame of Reference UID	(0020,0200)	1	A value generated by NW500.
Synchronization Trigger	(0018,106A)	1	NO TRIGGER
Acquisition Time Synchronized	(0018,1800)	1	N

A.2.6 General Equipment Module

Attribute Name	Tag	Type	Attribute Description
Manufacturer	(0008,0070)	2	A preset value
Institution Name	(0008,0080)	3	A preset value
Station Name	(0008,1010)	3	A preset value
Institutional Department Name	(0008,1040)	3	A preset value
Manufacturer's Model Name	(0008,1090)	3	NW500
Device Serial Number	(0018,1000)	3	Device Serial Number
Software Versions	(0018,1020)	3	Software Version

A.2.7 General Image Module

Attribute Name	Tag	Type	Attribute Description
Instance Number	(0020,0013)	2	A value generated by NW500.
Patient Orientation	(0020,0020)	2C	L¥F
Content Date	(0008,0023)	2C	A value generated by NW500.
Content Time	(0008,0033)	2C	A value generated by NW500.
Image Type	(0008,0008)	3	ORIGINAL¥PRIMARY
Acquisition Number	(0020,0012)	3	A value generated by NW500. In case of OP, this tag exists in an IOD.
Acquisition Date	(0008,0022)	3	Same as Content Date.
Acquisition Time	(0008,0032)	3	Same as Content Time.
Acquisition Date Time	(0008,002A)	3	A value generated from Content Date and Content Time. In case of OP, this tag exists in an IOD.
Burned In Annotation	(0028,0301)	3	NO
Lossy Image Compression	(0028,2110)	3	00 or 01

A.2.8 Image Pixel Module

Attribute Name	Tag	Type	Attribute Description
Sample per Pixel	(0028,0002)	1	A value generated by NW500.
Photometric Interpretation	(0028,0004)	1	A value generated by NW500.
Rows	(0028,0010)	1	A value generated by NW500.
Columns	(0028,0011)	1	A value generated by NW500.
Bits Allocated	(0028,0100)	1	A value generated by NW500.
Bits Stored	(0028,0101)	1	A value generated by NW500.
High Bit	(0028,0102)	1	A value generated by NW500.
Pixel Representation	(0028,0103)	1	A value generated by NW500.
Pixel Data	(7FE0,0010)	1C	A value generated by NW500.
Planar Configuration	(0028,0006)	1C	A value generated by NW500.

A.2.9 Acquisition Context Module

Attribute Name	Tag	Type	Attribute Description
Acquisition Context Sequence	(0040,0555)	2	Empty. In case of OP and VL, this tag is used.

A.2.10 SOP Common Module

Attribute Name	Tag	Type	Attribute Description
SOP Class UID	(0008,0016)	1	A value generated by NW500.
SOP Instance UID	(0008,0018)	1	A value generated by NW500.
Specific Character Set	(0008,0005)	1C	In case of using Japanese, this tag is used.
Instance Creation Date	(0008,0012)	3	A value generated by NW500.
Instance Creation Time	(0008,0013)	3	A value generated by NW500.

A.2.11 SC Equipment Module

In case of 'Secondary Capture Image', this module is used.

Attribute Name	Tag	Type	Attribute Description
Conversion Type	(0008,0064)	1	WSD

A.2.12 VL Image Module

In case of 'VL Photographic Image', this module is used.

Attribute Name	Tag	Type	Attribute Description
Image Type	(0008,0008)	1	ORIGINAL \backslash PRIMARY
Photometric Interpretation	(0028,0004)	1	A value generated by NW500.
Bits Allocated	(0028,0100)	1	A value generated by NW500.
Bits Stored	(0028,0101)	1	A value generated by NW500.
High Bit	(0028,0102)	1	A value generated by NW500.
Pixel Representation	(0028,0103)	1	A value generated by NW500.
Sample per Pixel	(0028,0002)	1	A value generated by NW500.
Planar Configuration	(0028,0006)	1C	A value generated by NW500.
Content Time	(0008,0033)	1C	A value generated by NW500.
Lossy Image Compression	(0028,2110)	2	00 or 01
Anatomic Region Sequence	(0008,2218)	1C	-
> Code Value	(0008,0100)	2	5665001
> Coding Scheme Designator	(0008,0102)	2	SCT
> Code Meaning	(0008,0104)	2	Retina
Pixel Data	(7FE0,0010)	1	A value generated by NW500.
Pixel Spacing	(0028,0030)	3	1.85E-03 \backslash 1.85E-03

A.2.13 SC Image Module

In case of 'Secondary Capture Image', this module is used.

Attribute Name	Tag	Type	Attribute Description
Pixel Spacing	(0028,0030)	1C	1.85E-03¥1.85E-03

A.2.14 Cine Module

In case of 'Ophthalmic Photography 8 Bit Image', this module is used.

Attribute Name	Tag	Type	Attribute Description
Frame Time	(0018,1063)	1C	0

A.2.15 Multi Frame Module

In case of 'Ophthalmic Photography 8 Bit Image', this module is used.

Attribute Name	Tag	Type	Attribute Description
Number of Frames	(0028,0008)	1	1
Frame Increment Pointer	(0028,0009)	1	(0018,1063)

A.2.16 Ophthalmic Photography Image Module

In case of 'Ophthalmic Photography 8 Bit Image', this module is used.

Attribute Name	Tag	Type	Attribute Description
Image Type	(0008,0008)	1	ORIGINAL¥PRIMARY¥¥COLOR
Instance Number	(0020,0013)	1	A value generated by NW500.
Samples per Pixel	(0028,0002)	1	A value generated by NW500.
Photometric Interpretation	(0028,0004)	1	A value generated by NW500.
Pixel Representation	(0028,0103)	1	A value generated by NW500.
Planar Configuration	(0028,0006)	1C	A value generated by NW500.
Pixel Spacing	(0028,0030)	1C	1.85E-03¥1.85E-03
Content Time	(0008,0033)	1	A value generated by NW500.
Content Date	(0008,0023)	1	A value generated by NW500.
Acquisition Date Time	(0008,002A)	1C	A value generated by NW500.
Lossy Image Compression	(0028,2110)	1	00 or 01
Burned In Annotation	(0028,0301)	1	NO

A.2.17 Ocular Region Imaged Module

In case of 'Ophthalmic Photography 8 Bit Image', this module is used.

Attribute Name	Tag	Type	Attribute Description
Image Laterality	(0020,0062)	1	R, L or B.
Anatomic Region Sequence	(0008,2218)	1	-
> Code Value	(0008,0100)	2	5665001
> Coding Scheme Designator	(0008,0102)	2	SCT
> Code Meaning	(0008,0104)	2	Retina

A.2.18 Ophthalmic Photography Acquisition Parameters Module

In case of 'Ophthalmic Photography 8 Bit Image', this module is used.

Attribute Name	Tag	Type	Attribute Description
Patient Eye Movement Commanded	(0022,0005)	2	YES or NO
Patient Eye Movement Commanded Code Sequence	(0022,0006)	1C	If (0022,0005) tag has YES, this tag is used.
> Code Value	(0008,0100)	2	A value defined in CID 4201 Patient Eye Movement command.
> Coding Scheme Designator	(0008,0102)	2	A value defined in CID 4201 Patient Eye Movement command.
> Code Meaning	(0008,0104)	2	A value defined in CID 4201 Patient Eye Movement command.
Refractive State Sequence	(0022,001B)	2	Empty
Emmetropic Magnification	(0022,000A)	2	A value generated by NW500.
Intra Ocular Pressure	(0022,000B)	2	Empty
Pupil Dilated	(0022,000D)	2	Empty

CID 4201 Patient Eye Movement Command

Coding Scheme Designator	Code Value	Code Meaning
SCT	408744005	Primary gaze
SCT	255533007	Upward gaze
SCT	255525006	Left upgaze
SCT	255530005	Left gaze
SCT	255523004	Left downgaze
SCT	255521002	Downgaze
SCT	255524005	Right downgaze
SCT	255531009	Right gaze
SCT	255526007	Right upgaze
SCT	408745006	Convergent gaze

A.2.19 Ophthalmic Photographic Parameters Module

In case of 'Ophthalmic Photography 8 Bit Image', this module is used.

Attribute Name	Tag	Type	Attribute Description
Acquisition Device Type Code Sequence	(0022,0015)	1	-
> Code Value	(0008,0100)	2	409898007
> Coding Scheme Designator	(0008,0102)	2	SCT
> Code Meaning	(0008,0104)	2	Fundus Camera
Detector Type	(0022,0016)	2	Empty
Light Path Filter Type Stack Code Sequence	(0022,0017)	2	Empty
Image Path Filter Type Stack Code Sequence	(0022,0018)	2	Empty
Lenses Code Sequence	(0022,0019)	2	Empty
Detector Type	(0018,7004)	2	Empty

End of report.