



OS-100



OS-200



OS-300



OS-400

Four different stool styles
to meet every need.

KEY FEATURES:


- » Extended pneumatic elevation (8 inches)
- » Premium casters that roll well on carpet and tile
- » Upholstered with bacteriostatic fabric
- » Sturdy built for comfort and durability
- » Easy to reach height adjustment lever

Performance & Specifications

Dimensions	OS-100	OS-200	OS-300	OS-400
Seat Size	15in (381mm)	15in (381mm)	16in x 15in (406mm x 381mm)	18in x 17in (457mm x 432mm)
Seat Foam Thickness	3in (76mm)	3in (76mm)	4in (102mm)	3in (76mm)
Retracted Seat Height	19in (483mm)	19in (483mm)	20in (507mm)	19in (483mm)
Extended Seat Height	27in (686mm)	27in (686mm)	28in (712mm)	27in (686mm)
Back Size	NA	10in x 15in (254mm x 381mm)	10in x 15in (254mm x 381mm)	12in x 17in (305mm x 432mm)
Base Footprint	24in x 24in (610mm x 610mm)	24in x 24in (610mm x 610mm)	24in x 24in (610mm x 610mm)	24in x 24in (610mm x 610mm)
Elevation Control	Side lever	Side lever	Side lever	Side lever
Caster Type	2in dual	2in dual	2.5in dual	2in dual/locking
Weight Capacity	250lbs (113.4 kg)	250lbs (113.4 kg)	250lbs (113.4 kg)	250lbs (113.4 kg)
Shipping Weight	21lbs (9.5kg)	28lbs (12.7kg)	31lbs (14kg)	31lbs (14kg)
Movements				
Elevation Stroke	8in (203mm)	8in (203mm)	8in (203mm)	8in (203mm)
Back Adjustment Range	NA	3in (76mm)	3in (76mm)	3in (76mm)
Appearance				
Upholstery Color	Your choice	Your choice	Your choice	Your choice
Cylinder Color	Black	Black	Chrome	Chrome
Base Color	Black	Black	Chrome	Chrome
Caster Color	Black	Black	Chrome	Chrome

Upholstery Selections

STANDARD COLORS



DESIGNER SERIES COLORS



Topcon Medical Systems, Inc.

111 Bauer Drive, Oakland, NJ 07436
 Phone: 800-223-1130 • Fax: 201-599-5248
 topconmedical.com

©2016 Topcon Medical Systems, Inc.



YOUR VISION. OUR FOCUS.

